

CSMHuman00000000033

Ultrasonic Doppler Vibrometry and the Spectral Signatures of Hemodynamic Turbulence: A Comprehensive Analysis of Frequency-Domain Human Body Interactions in Coronary Pathology

1. Introduction: The Paradigm Shift from Luminology to Hemodynamic Function

The detection and characterization of coronary artery disease (CAD) have traditionally been dominated by anatomical imaging modalities—a practice often pejoratively termed "luminology." Techniques such as X-ray angiography and intravascular ultrasound (IVUS) excel at visualizing the vessel lumen and the structural encroachment of atherosclerotic plaque. However, the geometric severity of a stenosis is an imperfect surrogate for its physiological impact. A lesion that appears anatomically moderate may induce severe hemodynamic compromise, while a seemingly severe calcification might be bypassed by robust collateral circulation, rendering it functionally benign. This discordance has driven the development of functional metrics like Fractional Flow Reserve (FFR), yet these typically require invasive catheterization.

This report explores a transformative non-invasive modality: Ultrasonic Doppler Vibrometry (UDV). By pivoting from the visualization of anatomy to the quantification of mechanical vibration, UDV interrogates the intrinsic fluid dynamics of the cardiovascular system.¹ The fundamental premise is that blood flow through a stenosis transitions from a laminar to a turbulent state, shedding vortices that impact the vessel wall. These impacts generate broad-spectrum mechanical vibrations—shear waves—that propagate through the viscoelastic myocardium and can be detected by high-frame-rate ultrasound.

This analysis synthesizes the physics of turbulent flow, the rheology of shear wave propagation in biological tissues, and the advanced signal processing architectures required to isolate these subtle signals from the high-amplitude noise of the beating heart. It further extends into the broader context of frequency-domain human body interactions, examining how these principles apply to passive elastography and pulmonary acoustics, ultimately proposing a unified theory of "spectral diagnostics" for soft tissue pathology.

2. Fluid Dynamics of Arterial Stenosis: The Genesis of Vibration

To understand the signal detected by UDV, one must first rigorously characterize the source: the hydrodynamic instability of poststenotic blood flow. The generation of a "bruit" or vibration is not a random artifact but a deterministic consequence of energy dissipation in fluid mechanics.

2.1 The Laminar-Turbulent Transition

In healthy coronary arteries, blood flow is predominantly laminar. The velocity profile is parabolic (Poiseuille flow), with maximal velocity at the center and zero velocity at the wall due to the no-slip condition. Viscous forces, governed by blood viscosity (μ), dampen any inertial perturbations. The stability of this flow is described by the Reynolds number (Re), a dimensionless ratio of inertial forces to viscous forces:

$$Re = \frac{\rho v D}{\mu}$$

Where ρ is blood density, v is mean velocity, and D is the vessel diameter. In normal coronary arteries, Re typically ranges from 200 to 500, well below the critical threshold for turbulence ($Re_{crit} \approx 2000-2300$ for pipe flow).³

However, the introduction of a stenosis radically alters this equilibrium. As blood enters the narrowed throat of a lesion, convective acceleration occurs to satisfy the conservation of mass ($A_1 v_1 = A_2 v_2$). This acceleration stabilizes the flow within the throat itself. The instability arises *distal* to the stenosis, in the "poststenotic dilation" or normal vessel segment. Here, the high-velocity jet exits into a wider lumen, creating an adverse pressure gradient. The boundary layer separates from the wall, creating a free shear layer between the high-velocity central jet and the stagnant or recirculating fluid near the walls.³

Research indicates that even in mild-to-moderate stenosis (30-50%), the flow remains largely laminar or transitional. However, as the stenosis severity exceeds 50%, the Reynolds number in the jet can exceed critical values, leading to the breakdown of the shear layer into coherent vortex structures.⁴ This transition is not binary but spectral: it evolves from stable laminar flow, to laminar flow with periodic vortex shedding, to fully developed broadband turbulence.

2.2 Vortex Shedding and the Strouhal Number

The mechanism most relevant to UDV is vortex shedding. The shear layer rolls up into vortices that detach periodically and convect downstream. When these vortices impact the arterial wall or interact with the slow-moving fluid, they generate pressure fluctuations. These

fluctuations are the acoustic source of the coronary bruit.

The frequency (f) of this vortex shedding is governed by the Strouhal number (St), a dimensionless parameter describing oscillating flow mechanisms:

$$St = \frac{f D}{v}$$

For a simple cylinder in cross-flow, St is approximately 0.2 across a wide range of Reynolds numbers. However, the complex geometry of an arterial stenosis—often eccentric and biologically compliant—alters this relationship. Extensive in vitro and in silico studies suggest that for arterial stenosis, the dominant "break frequency" corresponds to a Strouhal number significantly higher, often cited around 0.93.⁶

This finding is critical for sensor design. Given typical coronary flow velocities in a stenosis ($v \approx 100 - 400$ cm/s) and vessel diameters ($D \approx 0.1 - 0.3$ cm), the shedding frequency f falls into the range of:

$$f = \frac{0.93 \cdot 200 \text{ cm/s}}{0.2 \text{ cm}} \approx 930 \text{ Hz}$$

This calculation aligns with empirical observations that the energy of coronary bruits is concentrated in the 100 Hz to 1000 Hz band.¹ This frequency band is distinct from the fundamental heart rate (~1-2 Hz) and the wall motion associated with contraction (~20-30 Hz), allowing for spectral separation. However, it requires ultrasound systems capable of sampling frequencies well above the Nyquist limit of standard color Doppler (which typically operates at PRFs of 2-4 kHz, theoretically sufficient but practically limited by clutter filters).

2.3 The "Inverted U" Energy Paradox

One of the most nuanced and clinically significant findings in UDV research is the non-linear relationship between stenosis severity and vibration energy. Intuitively, one might expect vibration energy to increase monotonically with stenosis severity. However, experimental data from Sikdar, Beach, and Kim reveal a biphasic or "inverted U" relationship.¹

- **Phase 1: Acceleration (0% to ~60% Stenosis):** As the lumen narrows, jet velocity (v) increases. Since kinetic energy scales with v^2 and turbulent production scales with velocity gradients, the vibration energy rises sharply. The flow rate (Q) is maintained by distal autoregulation (dilation of the microvasculature).
- **Phase 2: The Peak (~60% to 70% Stenosis):** The vibration energy reaches its maximum. The velocity is high, and the flow rate is still sufficient to feed the turbulent structures.
- **Phase 3: Attenuation (>70% to 99% Stenosis):** As the stenosis becomes critical, the resistance (R) rises asymptotically ($R \propto 1/r^4$). Despite maximal vasodilation downstream, the pressure drop is so severe that the total volume flow rate (Q) begins to plummet.

The turbulent kinetic energy (\$TKE\$) generated is roughly proportional to the product of the flow rate and the pressure drop:

$Power_{turb} \propto Q \cdot \Delta P$

In critical stenosis, the decline in \$Q\$ outweighs the increase in \$\Delta P\$. The "mass" of blood available to participate in the turbulence decreases. Consequently, a critically occluded artery may become "quiet" or exhibit significantly reduced vibration energy compared to a moderate stenosis.⁸ This thermodynamic reality creates a potential diagnostic blind spot (the "quiet" vessel could be normal or totally occluded), necessitating the integration of UDV with anatomical imaging to resolve the ambiguity.

Table 1: Hemodynamic Regimes and Spectral Signatures

Stenosis Severity	Reynolds Number (Re)	Flow Regime	Flow Rate (Q)	Pressure Drop (ΔP)	Vibration Energy (DVE)
Normal (<30%)	< 500	Laminar / Poiseuille	Normal	Negligible	Baseline / Noise Floor
Mild (30-50%)	500 - 1500	Disturbed Laminar	Normal	Low	Slight Elevation
Moderate (50-70%)	1500 - 2500	Transitional / Vortex Shedding	Normal (Compensated)	Moderate	Maximum (Peak Signal)
Severe (70-90%)	> 2500	Fully Turbulent	Reduced	High	Decreasing
Critical (>95%)	Indeterminate	Trickle Flow / Relaminarization	Severely Reduced	Maximal	Low / Absent

3. Aeroacoustics and Fluid-Structure Interaction (FSI)

The turbulence described above occurs within the fluid. For UDV to detect it, the energy must

transfer from the blood to the vessel wall and then propagate through the surrounding tissue. This involves the physics of fluid-structure interaction (FSI) and aeroacoustics.

3.1 The Quadrupole Source and Lighthill's Analogy

Acoustically, turbulence in a free field is modeled as a quadrupole source, as defined by Lighthill's acoustic analogy. Quadrupole sources are notoriously inefficient radiators of sound, with acoustic power scaling as the eighth power of velocity (v^8).¹¹ However, in the cardiovascular system, the turbulence is not in a free field; it is bounded by a compliant vessel wall.

The interaction between the turbulent pressure fluctuations and the solid boundary enhances the radiation efficiency, effectively converting the source to a dipole character (v^6 scaling) or even causing direct mechanical coupling.⁶ The fluctuating pressure field acts as a forcing function on the arterial wall, exciting vibration modes. The vessel wall does not merely reflect the sound; it acts as an aeroacoustic transducer, converting hydrodynamic energy into mechanical displacement waves.

3.2 Shear Wave Generation

The vibration of the arterial wall generates two primary types of waves in the surrounding myocardium: compression waves (P-waves) and shear waves (S-waves).

- **Compression Waves:** These are longitudinal waves where particle motion is parallel to the direction of propagation. In soft tissue, they travel at approximately 1540 m/s. They are efficient carriers of energy over long distances but offer poor localization because they energize the entire thorax.
- **Shear Waves:** These are transverse waves where particle motion is perpendicular to propagation. In the myocardium, they travel much slower, typically 1 to 10 m/s, depending on stiffness.¹²

The distinction is crucial for UDV. The turbulence primarily excites *shear* deformation in the vessel wall due to the viscous drag and wall-normal pressure fluctuations. These shear waves propagate outward into the perivascular tissue. Unlike compression waves, shear waves in soft tissue are highly dispersive and heavily attenuated at high frequencies.

3.3 Frequency-Dependent Attenuation and Localization

The attenuation of shear waves in biological tissue follows a power law, often modeled as:

$$\alpha(f) = \alpha_0 f^n$$

For shear waves in the frequency band of coronary bruits (100–1000 Hz), the attenuation is extreme—often exceeding 10–20 dB per centimeter.¹⁴

This high attenuation, while seemingly a disadvantage, is actually the physical basis for the spatial specificity of UDV. A 500 Hz shear wave generated by the Left Anterior Descending (LAD) artery will decay to noise levels within 1-2 centimeters of the vessel. Therefore, if a focused ultrasound beam detects a 500 Hz vibration signal at a specific depth, it *must* originate from a source in the immediate vicinity (i.e., the coronary artery). It cannot be "cross-talk" from a valve or the aorta located 5 cm away, as those shear waves would have dissipated. This "proximity effect" allows UDV to map stenosis location with high spatial resolution without complex beamforming arrays required for source localization in lower-frequency seismocardiography.¹⁶

4. Wave Propagation in Viscoelastic Myocardium

The medium through which these diagnostic vibrations travel—the myocardium—is a complex material. It is anisotropic, viscoelastic, and actively changing its properties throughout the cardiac cycle. Understanding the rheology of the myocardium is essential for interpreting the UDV signal and links this technology to the broader field of elastography.

4.1 Viscoelastic Models: Kelvin-Voigt and Zener

Biological tissues exhibit both elastic (energy storage) and viscous (energy dissipation) behavior. The propagation of shear waves is mathematically described by the complex shear modulus (G^*):

$$G^*(\omega) = G'(\omega) + iG''(\omega)$$

Where G' is the storage modulus (elasticity) and G'' is the loss modulus (viscosity). The simplest model applied in UDV research is the Kelvin-Voigt model¹², which represents the tissue as a spring and dashpot in parallel. The stress (σ) is related to strain (ϵ) by:

$$\sigma = \mu_1 \epsilon + \mu_2 \frac{\partial \epsilon}{\partial t}$$

Here, μ_1 represents the static shear elasticity and μ_2 represents the shear viscosity. The phase velocity (c_s) of a shear wave in a Kelvin-Voigt material is dispersive (frequency-dependent):

$$c_s(\omega) \approx \sqrt{\frac{2(\mu_1^2 + \omega^2 \mu_2^2)}{\rho(\mu_1 + \sqrt{\mu_1^2 + \omega^2 \mu_2^2})}}$$

At low frequencies, the speed is dominated by elasticity ($\sqrt{\mu_1/\rho}$). At the high frequencies utilized in UDV (>100 Hz), the viscous term ($\omega \mu_2$) becomes significant, increasing the apparent stiffness and wave speed. This dispersion must be

accounted for when attempting to correlate vibration frequency with tissue properties.¹³

4.2 Myocardial Anisotropy and Fiber Orientation

The myocardium is not isotropic; it consists of muscle fibers arranged in a helical structure that varies transmurally from the endocardium to the epicardium. Shear waves propagate faster along the fiber direction than across it. Research using high-frequency shear wave elastography has shown that wave speeds can vary by a factor of 2 depending on the angle relative to the fiber orientation.¹³

For UDV, this implies that the detected vibration amplitude depends not only on the stenosis severity but also on the orientation of the ultrasound beam relative to the myocardial fibers adjacent to the artery. A beam perpendicular to the fibers might detect a different spectral signature than one parallel to them. Advanced modeling, potentially involving transversely isotropic viscoelastic layers, is required to fully deconvolve these effects.¹³

4.3 Passive Elastography: The Body as a Seismometer

The principles of UDV overlap significantly with "Passive Elastography." In active elastography (ARFI), an external force is applied. In passive elastography, the "noise" of the body—heartbeats, muscle twitches, and indeed, vascular turbulence—serves as the excitation source.

By analyzing the time-delay and correlation of shear waves generated by these internal sources, one can reconstruct the stiffness map of the tissue. This concept, borrowed from seismology (where ambient ground noise is used to map geological strata), suggests that UDV data contains latent information about myocardial stiffness. A stiff, fibrotic myocardium (common in ischemic heart disease) would support faster shear wave propagation with different attenuation characteristics than healthy tissue.¹⁸ Thus, UDV simultaneously interrogates the vessel (source) and the tissue (medium).

5. Ultrasonic Doppler Vibrometry (UDV): Instrumentation and Acquisition

Translating the physics of shear waves into a clinical image requires specialized ultrasonic instrumentation. Standard clinical scanners are optimized for B-mode (amplitude detection) and Color Doppler (mean frequency shift of blood). UDV requires access to the raw phase information at high temporal resolution.

5.1 The Doppler Principle Applied to Vibration

In UDV, the target is not the flowing blood, but the vibrating tissue. When an ultrasound pulse

with center frequency f_0 reflects off a vibrating target with displacement $x(t)$, the received signal's phase $\phi(t)$ is modulated:

$$\phi(t) = \frac{4\pi f_0}{c} x(t)$$

$$x(t) = A_{\text{vib}} \sin(2\pi f_{\text{vib}} t)$$

The velocity of the tissue $v(t)$ is the time derivative of displacement. The resulting Doppler shift f_d is:

$$f_d(t) = \frac{2 v(t) f_0}{c}$$

Since the vibration amplitude A_{vib} is microscopic (on the order of 0.1 to 10 μm) and the frequency f_{vib} is high (up to 1000 Hz), the resulting phase shifts are minute.

Detecting them requires a system with exceptionally low phase noise and high stability.¹

5.2 The Ultrasound Research Interface (URI)

Most clinical validations of UDV have utilized scanners equipped with a research interface, such as the Philips URI or the Verasonics vantage platform. These interfaces allow the user to bypass the standard backend processing (scan conversion, video compression) and access the raw In-Phase and Quadrature (IQ) data.²⁰

Key Acquisition Requirements:

1. **Pulse Repetition Frequency (PRF):** To detect vibrations up to 1000 Hz, the Nyquist theorem dictates a sampling rate of at least 2000 Hz. However, to accurately resolve the waveform for time-frequency analysis, oversampling is preferred. Typical cardiac PRFs of 4-8 kHz are sufficient.
2. **Dwell Time:** Standard color Doppler uses a packet size of 8-16 pulses per line to estimate mean velocity. This is insufficient for spectral analysis. UDV requires a "long dwell" acquisition, where the beam is held stationary (M-mode like) or scanned very slowly, acquiring hundreds of milliseconds of data per location to capture the full diastolic interval.¹
3. **Beamforming:** The range gate is positioned on the perivascular tissue, not the lumen. This avoids the "blooming" artifact of blood flow and focuses on the mechanical response of the wall.

Table 2: Comparison of Standard Doppler vs. UDV

Parameter	Standard Color Doppler	Ultrasonic Doppler Vibrometry (UDV)

Target	Red Blood Cells	Vessel Wall / Perivascular Tissue
PRF	2 - 4 kHz	4 - 10 kHz
Packet Size	8 - 16 pulses	> 200 pulses (Continuous)
Output	Mean Velocity, Variance	Displacement Time Series $x(t)$
Wall Filter	High-Pass (removes tissue motion)	Band-Pass (isolates tissue vibration)
Clutter Management	Rejects high-amplitude signals	Extracts high-freq/low-amp components
Clinical Goal	Lumen Patency / Regurgitation	Turbulence Quantification

6. Advanced Signal Processing Architectures

The raw IQ data contains a mixture of signals: the high-amplitude, low-frequency motion of the beating heart (clutter), the backscatter from blood, and the low-amplitude, high-frequency vibration of interest. Separating these components is the central challenge of UDV signal processing.

6.1 Eigen-Based Clutter Filtering (SVD)

Standard high-pass filters (IIR/FIR) are often inadequate because the spectral overlap between the rapid diastolic expansion of the ventricle and the low-end of the turbulent bruit (around 100 Hz) is significant. Furthermore, the "clutter" (heart wall motion) is non-stationary.

The state-of-the-art solution involves Singular Value Decomposition (SVD). The spatiotemporal data matrix $\mathbf{S}$ is constructed from the IQ samples.

$$\mathbf{S} = \mathbf{U} \mathbf{\Sigma} \mathbf{V}^*$$

In this decomposition:

- High Singular Values:** Represent the highly correlated, high-energy motion of the

myocardial wall (global contraction/relaxation).

- **Low Singular Values:** Represent random noise.
- **Intermediate Singular Values:** Typically contain the blood flow and the independent vibration components.

By zeroing out the first few singular values (which capture >90% of the energy associated with wall motion) and reconstructing the matrix, one can effectively perform a "smart" high-pass filter that adapts to the specific motion characteristics of the tissue.²¹ This technique, pioneered in ultrafast ultrasound, preserves the subtle shear wave information that would be destroyed by conventional filters.

6.2 Time-Frequency Analysis: The Morlet Wavelet

The vibration signal is transient, occurring primarily during the phase of maximal coronary flow (early diastole). A standard Fast Fourier Transform (FFT) applied over the whole cardiac cycle would dilute this signal with systolic noise.

The Continuous Wavelet Transform (CWT) using the **Morlet Wavelet** is the preferred method for spectral quantification.¹ The Morlet wavelet, being a Gaussian-windowed sinusoid, offers optimal time-frequency localization. It allows the algorithm to detect "bursts" of vibration energy—short-lived events (<100 ms) that characterize the turbulent shedding process.

The metric derived is the **Diastolic Vibration Energy (DVE)**, calculated by integrating the wavelet power spectrum over the 100–1000 Hz band specifically during the diastolic window.²³

$$DVE = \int_{t_{start}}^{t_{end}} \int_{100Hz}^{1000Hz} |CWT(t, f)|^2 df dt$$

6.3 Parametric Spectral Estimation

Beyond wavelets, researchers have explored parametric methods like Autoregressive (AR) and Autoregressive Moving Average (ARMA) models. These methods model the signal as the output of a linear system driven by white noise. While computationally more intensive, they can provide higher spectral resolution than FFT for short data segments. However, Sikdar et al. found that **Eigenvector methods** (such as MUSIC or ESPRIT) often outperformed FFT and AR methods in diagnostic accuracy, likely due to their ability to separate signal subspaces from noise subspaces effectively.¹

7. Clinical Evidence and Diagnostic Performance

The theoretical and engineering frameworks described above have been validated in several prospective clinical trials, primarily driven by the University of Washington and associated

groups.

7.1 Sensitivity and Specificity Profiles

In a landmark study of 49 patients undergoing coronary angiography, UDV demonstrated a robust ability to detect CAD.

- **Sensitivity:** 91.7% for detecting stenosis >50%.
- **Specificity:** 83.3%.¹⁶
- **ROC Analysis:** The Area Under the Curve (AUC) was 0.84, indicating good diagnostic classification performance.¹

These values are comparable to stress echocardiography and nuclear perfusion imaging (SPECT), but with the distinct advantage of being performed at rest, without the need for exercise or pharmacological stress agents.

7.2 Validation of the "Inverted U"

The clinical data strongly supported the "Inverted U" hypothesis derived from fluid dynamics. Patients with minor stenosis (<50%) showed low vibration energy (laminar flow). Patients with moderate stenosis (50-70%) exhibited the highest vibration energies (peak turbulence). Crucially, patients with severe stenosis (>70-90%) showed a statistically significant *decrease* in vibration energy compared to the moderate group.¹

This finding has major clinical implications:

1. **Screening Utility:** The test is highly sensitive for the "vulnerable" population with moderate stenosis who may not yet be symptomatic but are at risk of plaque rupture.
2. **False Negatives:** A "negative" UDV result could theoretically imply a critical occlusion. However, in clinical practice, a critical occlusion would be evident on standard B-mode echo (wall motion akinesis) or ECG. The true value of UDV is in the "grey zone" of intermediate lesions where standard imaging is ambiguous.

7.3 Differential Diagnosis and Artifacts

The chest is a noisy acoustic environment. The lungs create air-tissue interfaces that block ultrasound, and the heart valves (mitral and aortic) generate their own mechanical vibrations (S1, S2 heart sounds).

- **Valve Clicks:** These are high-amplitude, broadband transients. They can be filtered out by time-gating (analyzing only mid-diastole) and by their lower fundamental frequency (<50 Hz).
- **Lung Noise:** Respiration creates broad motion artifacts. These are typically low frequency (<5 Hz) and are removed by the SVD clutter filter. However, lung pathology itself (like pneumothorax) changes the transmission of sound to the chest wall, an interaction detailed in parallel research on lung acoustics.²⁴

8. Broader Applications: Passive Elastography and Pulmonary Acoustics

The principles of detecting high-frequency body vibrations extend beyond the coronary arteries. The interaction of frequency and human tissue is a fertile ground for "spectral diagnostics."

8.1 Passive Elastography and Time Reversal

As noted in Section 4.3, the shear waves generated by the heart propagate throughout the thorax. Recent advances in Passive Elastography utilize these waves to image tissue stiffness without an active "push" pulse.

By using Time Reversal algorithms or noise correlation interferometry (techniques adapted from seismology), researchers can virtually focus these diffuse shear waves to a point. If the focal spot is in a stiff region (e.g., a tumor or fibrosis), the wave speed and wavelength will be distinct.¹⁸ This suggests that the continuous physiological vibration of the human body is not noise, but a rich data source waiting to be decoded.

8.2 Pulmonary Vibration Analysis

The transmission of sound through the lungs is another example of frequency-body interaction. The lung parenchyma is a porous foam-like structure. Pathologies like fibrosis or pneumothorax radically alter its acoustic impedance and wave propagation speed.

Research using scanning laser Doppler vibrometry (a non-contact optical analog to UDV) on the chest wall has shown that:

- **Pneumothorax:** Causes a frequency-dependent decrease in sound transmission due to the air gap decoupling the lung from the chest wall.
- **Fibrosis:** Alters the shear wave motion but has negligible effect on compression waves. This highlights the importance of analyzing the *type* of wave (shear vs. compression) to diagnose specific tissue changes.²⁴

8.3 Seismocardiography (SCG)

While UDV looks deep inside the body, SCG measures the vibrations that reach the surface using accelerometers. Recent work has bridged these fields using "Airborne Ultrasound Surface Motion Cameras," which detect chest wall micro-vibrations remotely.²⁵ While SCG lacks the depth resolution of UDV, it operates on the same fundamental physics: the heart as a mechanical oscillator driving a viscoelastic medium.

9. Integration with Computational Fluid Dynamics

(CFD)

The future of UDV lies in its integration with Computational Fluid Dynamics (CFD). Currently, UDV provides a scalar measurement (energy) at a location. CFD can provide a vector field prediction.

By segmenting a patient's CT angiogram, researchers can create a 3D mesh of the coronary arteries. Solving the Navier-Stokes equations (using turbulence models like **k-omega SST** or **Large Eddy Simulation**) allows for the prediction of:

- 1. **Separation Zones:** Where turbulence will initiate.
- 2. **Wall Shear Stress (WSS):** The force acting on the endothelium.
- 3. **Pressure Fluctuations:** The acoustic source term.⁶

This creates a powerful feedback loop. The CFD model can predict the "break frequency" and expected vibration amplitude. If the UDV measurement matches the prediction, the diagnosis is confirmed. If they diverge, it implies the geometric model (CT) may have missed a lesion, or the tissue properties (stiffness) are abnormal. This "Digital Twin" approach allows for patient-specific calibration of the Strouhal number and turbulence thresholds.²⁸

Table 3: Turbulence Models in Computational Hemodynamics

Model	Complexity	Computational Cost	Applicability to Stenosis
Laminar (Steady)	Low	Low	Only for mild stenosis (<30%)
RANS (\$k-\epsilon\$)	Moderate	Low-Moderate	Good for mean flow, poor for separation
RANS (\$k-\omega\$ SST)	Moderate	Moderate	Standard for near-wall turbulence
Large Eddy Simulation (LES)	High	High	Captures vortex shedding / spectral content
Direct Numerical	Extreme	Extreme	Research only; resolves all

Sim (DNS)			turbulent scales
-----------	--	--	------------------

10. Conclusion and Future Outlook

Ultrasonic Doppler Vibrometry represents a convergence of three sophisticated disciplines: fluid dynamics, solid mechanics, and signal processing. By treating the human body not merely as a structure to be visualized, but as a medium of wave propagation and mechanical interaction, UDV unlocks a functional dimension of diagnostics.

The evidence presented in this report confirms that coronary stenosis generates a spectral signature—a "cry" of the artery—that is quantifiable, deterministic, and diagnostically potent. The transition to turbulence, governed by the Reynolds and Strouhal numbers, creates shear waves that broadcast the vessel's hemodynamic status through the myocardium.

While challenges remain—specifically the "Inverted U" ambiguity and the need for advanced clutter filtering—the trajectory of ultrasound technology favors UDV. The shift toward software-defined, ultrafast ultrasound systems (capable of >10,000 frames per second) removes the hardware bottlenecks that previously limited this technique.

Ultimately, UDV validates the concept that the frequency domain holds vital biological information. From the shear waves of coronary turbulence to the transmission of sound in fibrotic lungs, the analysis of human body interaction with mechanical frequency offers a non-invasive window into pathology that complements, and in some cases surpasses, the static images of luminology.

Report Authored By:
Domain Expert in Biomedical Engineering & Computational Hemodynamics
Specialization: Ultrasound Physics, Aeroacoustics, and Non-Invasive Diagnostics

Works cited

1. Ultrasonic Doppler Vibrometry: Novel Method for Detection of Left Ventricular Wall Vibrations Caused by Poststenotic Coronary Flow | Request PDF - ResearchGate, accessed December 23, 2025, https://www.researchgate.net/publication/6075324_Ultrasonic_Doppler_Vibrometry_Novel_Method_for_Detection_of_Left_Ventricular_Wall_Vibrations_Caused_by_Poststenotic_Coronary_Flow
2. Doppler vibrometry: assessment of arterial stenosis by using perivascular tissue vibrations without lumen visualization. - Semantic Scholar, accessed December 23, 2025, <https://www.semanticscholar.org/paper/Doppler-vibrometry%3A-assessment-of->

- [arterial-stenosis-Sikdar-Vaidya/5be4a6a75ec14efb290d6cf5464f65bec92c0294](#)
3. Detection of Moving Flow Separation in Pulsatile Flow and the Degree of Stenosis by Power of Doppler Shift Signals - American Heart Association Journals, accessed December 23, 2025, <https://www.ahajournals.org/doi/pdf/10.1161/01.res.67.1.166>
 4. Numerical analysis of the effect of turbulence transition on the hemodynamic parameters in human coronary arteries - ResearchGate, accessed December 23, 2025, https://www.researchgate.net/publication/341591685_Numerical_analysis_of_the_effect_of_turbulence_transition_on_the_hemodynamic_parameters_in_human_coronary_arteries
 5. Numerical analysis of the effect of turbulence transition on the hemodynamic parameters in human coronary arteries - PubMed Central, accessed December 23, 2025, <https://pmc.ncbi.nlm.nih.gov/articles/PMC4880750/>
 6. Computational Modeling and Analysis of Murmurs Generated by Modeled Aortic Stenoses - Johns Hopkins Whiting School of Engineering, accessed December 23, 2025, https://engineering.jhu.edu/fsag/wp-content/uploads/2019/05/JBME_2019.pdf
 7. An FSEI approach for the assessment of stenotic aortic valve effects on the left heart hemodynamics | Request PDF - ResearchGate, accessed December 23, 2025, https://www.researchgate.net/publication/372822701_An_FSEI_approach_for_the_assessment_of_stenotic_aortic_valve_effects_on_the_left_heart_hemodynamics
 8. Ultrasonic Doppler vibrometry: novel method for detection of left ventricular wall vibrations caused by poststenotic coronary flow - PubMed, accessed December 23, 2025, <https://pubmed.ncbi.nlm.nih.gov/17764895/>
 9. QT interval dispersion analysis in patients with aortic valve stenosis: A prospective study, accessed December 23, 2025, https://www.researchgate.net/publication/7004523_QT_interval_dispersion_analysis_in_patients_with_aortic_valve_stenosis_A_prospective_study
 10. Vascular Technology - Davies Publishing, accessed December 23, 2025, https://www.daviespublishing.com/assets/images/product_samples/sample-RPVI-Bund.pdf
 11. A CFD-FFT approach to hemoacoustics that enables degree of stenosis prediction from stethoscopic signals, accessed December 23, 2025, <https://pmc.ncbi.nlm.nih.gov/articles/PMC10336451/>
 12. VISCOELASTIC PROPERTIES OF NORMAL AND INFARCTED MYOCARDIUM MEASURED BY A MULTIFREQUENCY SHEAR WAVE METHOD: COMPARISON WITH PRESSURE-SEGMENT LENGTH METHOD - NIH, accessed December 23, 2025, <https://pmc.ncbi.nlm.nih.gov/articles/PMC4118646/>
 13. Modeling Transversely Isotropic, Viscoelastic, Incompressible Tissue-like Materials with Application in Ultrasound Shear Wave Elastography | Request PDF - ResearchGate, accessed December 23, 2025, https://www.researchgate.net/publication/270963528_Modeling_Transversely_Isotropic_Viscoelastic_Incompressible_Tissue-like_Materials_with_Application_in_Ultr

[asound_Shear_Wave_Elastography](#)

14. Longitudinal shear wave and transverse dilatational wave in solids - ResearchGate, accessed December 23, 2025, https://www.researchgate.net/publication/272386336_Longitudinal_shear_wave_and_transverse_dilatational_wave_in_solids
15. Attenuation Measuring Ultrasound Shearwave Elastography and in vivo application in post-transplant liver patients - PMC - NIH, accessed December 23, 2025, <https://pmc.ncbi.nlm.nih.gov/articles/PMC5505633/>
16. Ultrasonic Doppler Vibrometry: Measurement of Left Ventricular Wall Vibrations Associated With Coronary Artery Disease - IEEE Xplore, accessed December 23, 2025, <https://ieeexplore.ieee.org/iel5/4028925/4461641/04461888.pdf>
17. Dynamic Elastography of Pre-stressed Material with Multi-Directional Excitation in High Field MRI System - UIC Indigo, accessed December 23, 2025, https://indigo.uic.edu/articles/thesis/Dynamic_Elastography_of_Pre-stressed_Material_with_Multi-Directional_Excitation_in_High_Field_MRI_System/17026433/1/files/31491005.pdf
18. (PDF) Seismic time-lapse interferometry across scales - ResearchGate, accessed December 23, 2025, https://www.researchgate.net/publication/334643089_Seismic_time-lapse_interferometry_across_scales
19. 1-D Elasticity Assessment in Soft Solids from Shear Wave Correlation: The Time-Reversal Approach | Request PDF - ResearchGate, accessed December 23, 2025, https://www.researchgate.net/publication/40035681_1-D_Elasticity_Assessment_in_Soft_Solids_from_Shear_Wave_Correlation_The_Time-Reversal_Approach
20. Research interface on a programmable ultrasound scanner | Request PDF - ResearchGate, accessed December 23, 2025, https://www.researchgate.net/publication/5618187_Research_interface_on_a_programmable_ultrasound_scanner
21. Imaging Heart Dynamics With Ultrafast Cascaded-Wave Ultrasound - IEEE Xplore, accessed December 23, 2025, <https://ieeexplore.ieee.org/iel7/58/8813107/08746191.pdf>
22. Fast randomized singular value decomposition based clutter filtering for shear wave imaging | Request PDF - ResearchGate, accessed December 23, 2025, https://www.researchgate.net/publication/327880062_Fast_randomized_singular_value_decomposition_based_clutter_filtering_for_shear_wave_imaging
23. Ultrasonic Doppler vibrometry: measurement of left ventricular wall vibrations associated with coronary artery disease - PubMed, accessed December 23, 2025, <https://pubmed.ncbi.nlm.nih.gov/17946866/>
24. Sound transmission in human thorax through airway insonification: an experimental and computational study with diagnostic applications - PubMed Central, accessed December 23, 2025, <https://pmc.ncbi.nlm.nih.gov/articles/PMC7501255/>
25. Ondes et Images : Articles 2018 - Institut Langevin - ESPCI Paris, accessed December 23, 2025, https://www.institut-langevin.espci.fr/journal_articles_2018

26. Ondes et Images : Mathias FINK - Institut Langevin - ESPCI Paris, accessed December 23, 2025, https://www.institut-langevin.espci.fr/mathias_fink
27. Instantaneous Visual Analysis of Blood Flow in Stenoses Using Morphological Similarity, accessed December 23, 2025, <https://arxiv.org/html/2403.16653v1>
28. Effect of shape of the stenosis on the hemodynamics of a stenosed coronary artery - AIP Publishing, accessed December 23, 2025, https://pubs.aip.org/aip/pof/article-pdf/doi/10.1063/5.0058765/14132009/081914_1_online.pdf